







Sample changer overview

						
	862 Compact Titrator	869 Compact Sample Changer	814 USB Sample Processor	815 Robotic USB Sample Processor XL	855 Robotic Titrator	864 Robotic Balance Sample Processor
Integrated measuring input	yes	–	–	–	yes	–
Integrated balance	–	–	–	–	–	yes
Work stations	1	1	1 / 2	1 / 2	1	2
Automatic rack recognition	yes	yes	yes	yes	yes	yes
Exchangeable racks	–	–	yes	yes	yes	yes
Titration/measurement directly on the rack	yes	yes	yes	yes	yes	yes
Standard racks	11 x 120 mL	11 x 120 mL	12 x 250 mL 16 x 150 mL 22 x 120 mL 24 x 75 mL	28 x 250 mL 34 x 150 mL 59 x 120 mL 100 x 75 mL 160 x 6 mL	28 x 250 mL 34 x 150 mL 59 x 120 mL 100 x 75 mL 160 x 6 mL	20 x 120 mL
Thermostatic racks	–	–	–	yes	yes	–
Beaker sensor	–	–	yes	yes	yes	yes
Splash protection	yes	yes	yes	yes	yes	yes
Drip pan	–	–	yes	yes	yes	yes
Max. number internal pumps	–	–	2 / 4	2 / 4	2	–
Max. number external pumps	2	2	2 / 4	2 / 4	1	4
Connector Swing Head	–	–	1 / 2	1 / 2	1	2
External positions	–	–	4 / 8	4 / 8	4	8
Pipetting / Sample transfer	–	–	yes	yes	yes	yes
Stirrer connector	1	1	1 / 2	1 / 2	1	2
Connector MSB	1	–	3	3	3	3
Connector Remote	yes	yes	optional	optional	optional	optional
Connector USB	–	–	2	2	2	2
Connector Mini-USB	yes	yes	–	–	–	–
Stand-alone	yes	yes	optional	optional	optional	–
PC Control	–	–	yes	yes	yes	–
tiamo™	–	–	yes	yes	yes	yes