

Press Release: Dual Titrations using the Metrohm Titrando

Date: 17th August 2010

For the last year or so we have seen an increased trend in the demand for our new range of titrators, particularly titrations and dual titrations. This is not particular to a specific industry, whether the food industry, pharmaceutical chemical or contract analysis, we have supplied instruments that save time and space for titrations.

The Titrando range can perform a large range of applications and determinations, this article is not going to detail all of these but outline the dual titration concept.

The 905, 906 and 907 models can be used to perform dual and/or parallel titrations. What we mean by this is one titrator with two stirring stations performing independent titrations at the same time. This concept saves time, saves space as only one titrators is being used and best of all will save money.

The Titrando is a state of the art titrators and recent acquisitions for dual titrations have been made by the food industry with acidity and salt performed simultaneously, FFAs and peroxide values performed together; the plating industry with Cr III and Al content performed simultaneously and the pharmaceutical with aqueous and non-aqueous acid base determinations. Control couldn't be easier with the state-of-the-art Tiamo 2.2 controlling softwa



So if you want to improve your procedures, save time space and money contact us to arrange a demonstration on 01928 579 600 or email enquiry@metrohm.co.uk

For more information on the Metrohm Titrando range visit www.metrohm.co.uk

For more information on this press release please contact the Marketing Department at Metrohm UK Ltd.

Helen Davey - Marketing Manager

Metrohm UK Ltd

Direct Tel: 01928 532 386. Gen Tel. 01928 579600

Fax: 01928 532371

Mobile: 07773 023 766

Note to Editors:

Metrohm UK Ltd is the UK distributor of Metrohm AG, the leading manufacturer of Ion Analysis Instruments.

Metrohm is synonymous with potentiometric titration, Karl Fischer titration, ion chromatography and voltammetry solving your analytical chemistry problems with quality instrumentation and support services in the technologies of POTENTIOMETRIC TITRATION, KARL FISCHER TITRATION (moisture analysis), ION CHROMATOGRAPHY, VOLTAMMETRY (trace analysis), On-line & AT line analysis and VISCOMETRY. Set up in 1995, Metrohm UK serves the key markets of the UK including the chemical, petrochemical and pharmaceutical industries, as well as the automotive and food sectors.