

## 25 Years of **Metrohm IC**

Robust, Reliable, Responsive Ion Chromatography Systems

### **MiPT Upgrade Offer\*** (Metrohm intelligent Partial Loop Technique)



**5 YEAR  
INSTRUMENT  
WARRANTY**

**£2899\*\***

The Metrohm intelligent Partial Loop Technique (MiPT) uses our patented Dosino technology to deliver variable injection volumes to a Metrohm IC system.

\*MiPT can be added to any Metrohm IC system containing the 858 Professional Sample Processor.

#### Unique features:

- Automatic calibration using just one single multi-ion standard
- Injection volumes between 2 and 200  $\mu\text{L}$
- Smallest carryover from sample to sample of all sample preparation techniques:  $<0.001\%$
- Sample concentrations in the range 1:10'000 can be analysed reliably with automatic calibration
- Capable of following logical decisions which can automatically adjust the injection volume to the concentration of the sample
- Exact and reliable results, because the results are always within the calibration range
- Faster and more precise ( $n=10$ : RSD  $<0.1\%$ ) than manual dilution

#### Offer includes:

- MiPT equipment
- Onsite installation and training ( $\sim\frac{1}{2}$  day)
- An additional  $\frac{1}{2}$  day of the technical representative's time to be used as required e.g. additional staff basic training, advanced software training, troubleshooting, method optimisation etc.

For more details please contact us on 01928 579 600 or [info@metrohm.co.uk](mailto:info@metrohm.co.uk) where we will be delighted to provide further information upon request.

\*\*Offer available until the 31<sup>st</sup> March 2012  
5 year warranty only valid with proof of service history