

Press Release: Better handling: Metrohm Aquatrode plus now also with removable electrode cable

Date: 24th January 2012

The Aquatrode plus with plug-in head U from Metrohm was specially developed for measurements and titrations in low-ion and poorly buffered samples, such as surface waters and deionized water. These are the major benefits:



- Special membrane glass with low resistance for measurement in low-ion solutions. This guarantees very short response times even in such difficult samples (fast measurement or end point titration).
- Fixed ground-joint diaphragm with low resistance and very low, clearly defined electrolyte flow for fast, very accurate and reproducible pH measurement / titration.
- Double-junction design with maintenance-free inner electrolyte (gel) enables the bridge electrolyte to be adapted to the measurement conditions.
- Screening of the inner tube prevents interference due to static electricity
- In addition to these advantages, the Aquatrode plus now offers easy storage / fast electrode replacement / complete flexibility due to its removable cable.

For more information on this press release please contact the Marketing Department at Metrohm UK Ltd.

Helen Davey - Marketing Manager Metrohm UK Ltd
Direct Tel: 01928 532 386. Gen Tel. 01928 579600

Note to Editors:

Metrohm is synonymous with potentiometric titration, Karl Fischer titration, ion chromatography and electrochemistry solving your analytical chemistry problems with quality instrumentation and support services set up in 1995, Metrohm UK serves the key markets of the UK including the chemical, petrochemical and pharmaceutical industries, as well as the automotive and food sectors.