

Press Release: Space saving and robust: The new combined Ca-ISE

Date: 7th June 2011



No matter if titration, direct measurement or standard addition, the new combined ion-selective Calcium sensor is the ideal choice. It is especially suited for drinking water analysis where other parameters e.g. chloride have to be determined at the same time.

The combined sensor requires less space and comes with a chloride-free electrolyte – contamination of the sample with chloride is avoided.

To find out more visit www.metrohm.ie

For more information on this press release please contact the Marketing Department at Metrohm Ireland Ltd.

Helen Davey - Marketing Manager

Metrohm Ireland Ltd

Tel. 01928 579600

Email hdavey@metrohm.co.uk

Note to Editors:

Metrohm Ireland Ltd is the Irish distributor of Metrohm AG, the leading manufacturer of Ion Analysis Instruments.

Metrohm is synonymous with potentiometric titration, Karl Fischer titration, ion chromatography and voltammetry solving your analytical chemistry problems with quality instrumentation and support services in the technologies of POTENTIOMETRIC TITRATION, KARL FISCHER TITRATION (moisture analysis), ION CHROMATOGRAPHY, VOLTAMMETRY (trace analysis), On-line & AT line analysis and VISCOMETRY.

Set up in 1995, Metrohm UK serves the key markets of the UK including the chemical, petrochemical and pharmaceutical industries, as well as the automotive and food sectors.