

Potentiometric Titrators - overview

Potentiometrische Titratoren	Titrino Plus			Titrando							
	870	877	848	888	890	901	902	904	905	906	907
Optional 2nd measuring interface	-	-	-	-	-	yes	yes	yes	yes	yes	yes
Intelligent Sensor «iTrodes»	-	-	-	yes	-	yes	yes	yes	yes	yes	yes
Volumetric Karl Fischer Titration	yes	-	-	-	yes	yes	-	-	-	yes	yes
Ion Measurement / Standardaddition	-	-	-	-	-	-	-	yes	yes	yes	yes
pH Measurement	-	yes	yes	yes	-	yes	yes	yes	yes	yes	yes
Endpoint Titration	-	yes	yes	yes	-	yes	yes	yes	yes	yes	yes
Dynamic Equivalence Point Titration	-	-	yes	yes	-	-	-	yes	yes	yes	yes
Monotonic Equivalence Point Titration	-	yes	yes	yes	-	-	-	yes	yes	yes	yes
Liquid Handling	-	-	-	-	-	-	yes	yes	yes	yes	yes
STAT Titration	-	-	-	-	-	-	yes	-	-	yes	yes
Dosino Technology	-	-	-	-	-	yes	yes	-	yes	-	yes
Software Control	-	-	-	yes	yes	yes	yes	yes	yes	yes	yes
Touch Control	-	-	-	yes	yes	yes	yes	yes	yes	yes	yes
Live Display	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes
Report Printing / LIMS	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes
Balance Connection via RS232	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes	yes

