

Press Release: Free monograph «Voltammetric analysis methods in electroplating»

Date: 20th August 2010

Due to their high selectivity and low detection limits, voltammetric analysis methods are indispensable in surface finishing and PCB production. The application spectrum ranges from bath control to wastewater analysis. The present monograph is free of charge and provides a comprehensive overview of specific methods and also covers the determination of organic additives.

The first part of the monograph is dedicated to wastewater analysis in the galvanic industry using stripping voltammetric methods. Analytical methods for chromium, cyanide, nickel, cobalt, zinc, cadmium, lead, copper, tin and the complexing agents EDTA and NTA are described. Since wastewater almost always contains organic compounds, samples must be pre-processed prior to analysis. For this purpose, the authors recommend UV treatment using Metrohm's 705 UV Digester.



The second and more comprehensive part of the monograph addresses the voltammetric analyses of galvanic baths. The authors provide a systematic overview of the different bath types for each element as well as their primary and secondary components. The most important analytical methods are subsequently presented step-by-step. The practical value of the book is further enhanced by numerous references to other cost-free publications (application bulletins, application notes).

Get a comprehensive overview of voltammetric analysis methods in surface finishing by ordering your free personal copy from Metrohm UK.

For more information on this press release please contact the Marketing Department at Metrohm UK Ltd.

Helen Davey - Marketing Manager

Metrohm UK Ltd

Direct Tel: 01928 532 386. Gen Tel. 01928 579600

Fax: 01928 532371

Mobile: 07773 023 766

Note to Editors:

Metrohm UK Ltd is the UK distributor of Metrohm AG, the leading manufacturer of Ion Analysis Instruments.

Metrohm is synonymous with potentiometric titration, Karl Fischer titration, ion chromatography and voltammetry solving your analytical chemistry problems with quality instrumentation and support services in the technologies of POTENTIOMETRIC TITRATION, KARL FISCHER TITRATION (moisture analysis), ION CHROMATOGRAPHY, VOLTAMMETRY (trace analysis), On-line & AT line analysis and VISCOMETRY.

Set up in 1995, Metrohm UK serves the key markets of the UK including the chemical, petrochemical and pharmaceutical industries, as well as the automotive and food sectors.