

## Volumetric KF titrators – overview



Technical specifications	870 KF Titrino plus	890 Titrande	901 Titrande
Volumetric Karl Fischer titration (KFT)	yes	yes	yes
Titration to preset pH/mV endpoints (SET)	–	–	yes
pH calibration and pH measurements (CAL, MEAS)	–	–	yes
Real-time curve display	yes	yes	yes
Predefined methods	yes	yes	yes
Parallel titration	–	–	yes ( <i>tiamo</i> <sup>TM</sup> )
Dosing drive for titration	Exchange Unit	Exchange Unit	Dosing Unit
<b>Inputs, interfaces</b>			
Second, galvanically separated measuring interface	–	–	option
Automated reagent exchange with Dosino (option)	–	yes	yes
MSB-port for connection of further Dosinos (addition of auxiliary solution, ...)	–	3	4
MSB-port for connection of further Dosinos (Titration)	–	–	3
Attachment of stirrer	yes	yes	yes
Sample changer connection	–	yes	yes
Attachment of printer, balance, ...	yes, USB	yes, USB	yes, USB
Optional interfaces	optional RS 232C	yes, USB	yes, USB
<b>Software</b>			
<i>tiamo</i> <sup>TM</sup>	–	yes (1.3 or higher)	yes (2.0 or higher)
PC Control or Touch Control	–	yes	yes (6.0 or higher)